

XPO4 Rabbit Polyclonal Antibody



CAB15503

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

130kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

64328

Uniprot

Q9C0E2

Synonyms:

XPO4; exp4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human XPO4 (NP_071904.4).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
