

Product Information

Product SKU:	CAB15514	Gene ID:	79876	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the E1-like ubiquitin-activating enzyme family. This protein activates ubiquitin-fold modifier 1, a ubiquitin-like post-translational modifier protein, via the formation of a high-energy thioester bond. Alternative splicing results in multiple transcript variants. A pseudogene of this gene has been identified on chromosome 1.

Recommended Dilution:	WB,1:200 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	DEE44; EIEE44; SCAR24; THIFP1; UBE1DC1; UBA5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-404 of human UBA5 (NP_079094.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.