

THNSL1 Rabbit Polyclonal Antibody

CAB15515



Product Information

Product SKU:	CAB15515	Gene ID:	79896	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	83kDa	Conjugate:	Unconjugated
Calculated MW:	83kDa	Isotype:	IgG

Immunogen Information

Background: THNSL1 (Threonine Synthase Like 1) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include pyridoxal phosphate binding and threonine synthase activity. An important paralog of this gene is THNSL2.

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	TSH1; THNSL1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 534-743 of human THNSL1 (NP_079114.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.