

PARP9 Rabbit Polyclonal Antibody



CAB15526

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

80kDa/92kDa/96kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

83666

Uniprot

Q8IXQ6

Synonyms:

PARP9; ARTD9; BAL; BAL1; MGC:7868

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human PARP9 (NP_001139577.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
