

ZNF564 Rabbit Polyclonal Antibody



CAB15572

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

63kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

163050

Uniprot

Q8TBZ8

Synonyms:

ZNF564

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human ZNF564 (NP_659413.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
