

# MIA3 Rabbit Polyclonal Antibody



CAB15582

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

56kDa/87kDa/206kDa/213kDa

a

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:200 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

375056

#### Uniprot

Q5JRA6

#### Synonyms:

MIA3; ARNT; D320; TANGO; TANGO1; UNQ6077

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1450-1640 of human MIA3 (NP\_940953.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---