

ZBTB18 Rabbit Polyclonal Antibody



CAB15610

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

58kDa/59kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a C2H2-type zinc finger protein which acts a transcriptional repressor of genes involved in neuronal development. The encoded protein recognizes a specific sequence motif and recruits components of chromatin to target genes. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

10472

Uniprot

Q99592

Synonyms:

ZBTB18; C2H2-171; MRD22; RP58; TAZ-1; ZNF238

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 140-360 of human ZBTB18 (NP_006343.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
