

CAB15615

Product Information

Product SKU:	CAB15615	Gene ID:	4333680	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

**Background:** Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	RGB1; OsWD40-80; OsRGB1
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsRGB1
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.