

CAB15632

Product Information

Product SKU:	CAB15632	Gene ID:	1027	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	30kDa	Conjugate:	Unconjugated
Calculated MW:	22kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a cyclin-dependent kinase inhibitor, which shares a limited similarity with CDK inhibitor CDKN1A/p21. The encoded protein binds to and prevents the activation of cyclin E-CDK2 or cyclin D-CDK4 complexes, and thus controls the cell cycle progression at G1. The degradation of this protein, which is triggered by its CDK dependent phosphorylation and subsequent ubiquitination by SCF complexes, is required for the cellular transition from quiescence to the proliferative state. Mutations in this gene are associated with multiple endocrine neoplasia type IV (MEN4).

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:100 IF/ICC,1:50 - 1:200

Synonyms: KIP1; MEN4; CDKN4; MEN1B; P27KIP1; CDKN1B/p27KIP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDKN1B/p27KIP1 (NP_004055.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.