

CAB15692

Product Information

Product SKU:	CAB15692	Gene ID:	4664	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	54kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable transcription coregulator activity. Predicted to be involved in regulation of transcription, DNA-templated. Predicted to act upstream of or within endochondral ossification; nervous system development; and regulation of epidermis development. Predicted to be located in nucleoplasm. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NAB1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 189-317 of human NAB1 (NP_005957.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.