
Product Information

Product SKU:	CAB15725	Gene ID:	6617	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50-55kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable sequence-specific DNA binding activity. Predicted to be involved in snRNA transcription by RNA polymerase II and snRNA transcription by RNA polymerase III. Located in nucleolus and nucleoplasm.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: SNAP43; PTFgamma; SNAPC1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-368 of human SNAPC1 (NP_003073.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.