

Product Information

Product SKU:	CAB15738	Gene ID:	7692	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	73kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable DNA-binding transcription repressor activity, RNA polymerase II-specific and RNA polymerase II transcription regulatory region sequence-specific DNA binding activity. Involved in negative regulation of transcription by RNA polymerase II. Predicted to be located in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ZNF150; pHZ-13; pHZ-66; ZNF133
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ZNF133 (NP_003425.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.