

USP3 Rabbit Polyclonal Antibody

CAB15769



Product Information

Product SKU:	CAB15769	Gene ID:	9960	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	67kDa	Conjugate:	Unconjugated
Calculated MW:	59kDa	Isotype:	IgG

Immunogen Information

Background:	Enables histone binding activity and thiol-dependent deubiquitinase. Involved in several processes, including DNA repair; histone deubiquitination; and regulation of protein stability. Located in several cellular components, including Flemming body; cytoplasmic ribonucleoprotein granule; and nucleoplasm.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	UBP; SIH003; USP3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 321-520 of human USP3 (NP_006528.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.