

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB15775 | Gene ID: | 10200 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|-------------|------------|--------------|
| Observed MW: | 22kDa,23kDa | Conjugate: | Unconjugated |
| Calculated MW: | 19kDa | Isotype: | IgG |

Immunogen Information

Background: Predicted to enable RNA binding activity. Involved in maturation of 5.8S rRNA. Located in cytosol; nucleolus; and nucleoplasm. Colocalizes with nuclear exosome (RNase complex).

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Synonyms: | MPP; MPP6; MPP-6; MPHOSPH6 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human MPHOSPH6 (NP_005783.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |