

RABGAP1 Rabbit Polyclonal Antibody



CAB15804

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

122kDa

Calculated MW:

28kDa/53kDa/67kDa/121kDa

Applications:

WB IHC

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23637

Uniprot

Q9Y3P9

Synonyms:

RABGAP1; GAPCENA; TBC1D11


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 830-1069 of human RABGAP1 (NP_036329.3).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot analysis of extracts of various cell lines, using RABGAP1 antibody (CAB15804) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 90s.



Immunohistochemistry of paraffin-embedded Human colon carcinoma using RABGAP1 Rabbit pAb (CAB15804) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse testis using RABGAP1 Rabbit pAb (CAB15804) at dilution of 1:100 (40x lens).