

Product Information

Product SKU:	CAB15804	Gene ID:	23637	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	122kDa	Conjugate:	Unconjugated
Calculated MW:	122kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables GTPase activator activity and small GTPase binding activity. Involved in regulation of GTPase activity. Located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	GAPCENA; TBC1D11; RABGAP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 830-1069 of human RABGAP1 (NP_036329.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.