

Product Information

Product SKU:	CAB15807	Gene ID:	25940	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables protein methyltransferase activity. Involved in positive regulation of cell population proliferation; positive regulation of gene expression; and protein methylation. Predicted to be part of tRNA-splicing ligase complex.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:100 - 1:200

**Synonyms:** FAM98A

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human FAM98A (NP\_056290.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.