

Product Information

Product SKU:	CAB15809	Gene ID:	26009	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	170kDa	Conjugate:	Unconjugated
Calculated MW:	102kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable DNA binding activity and zinc ion binding activity. Predicted to be involved in histone H4 acetylation. Located in nucleolus and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ATAC1; ZZZ3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZZZ3 (Q8IYH5).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.