

CAB15828

Product Information

Product SKU:	CAB15828	Gene ID:	30849	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	150kDa	Conjugate:	Unconjugated
Calculated MW:	153kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable protein serine/threonine kinase activity. Involved in positive regulation of phosphatidylinositol 3-kinase activity; receptor catabolic process; and regulation of cytokinesis. Located in late endosome and microtubule cytoskeleton.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:100 IF/ICC,1:50 - 1:200

Synonyms: p150; VPS15; PIK3R4/VPS15

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human PIK3R4/VPS15 (NP_055417.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.