

PIK3R4 Rabbit Polyclonal Antibody



CAB15828

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

153kDa

Calculated MW:

153kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:100 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

30849

Uniprot

Q99570

Synonyms:

PIK3R4; VPS15; p150


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human PIK3R4 (NP_055417.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using PIK3R4 antibody (CAB15828) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 10s.



Immunofluorescence analysis of NIH-3T3 cells using PIK3R4 Polyclonal Antibody (CAB15828) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.