

TRIM62 Rabbit Polyclonal Antibody



CAB15855

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

54kDa

Calculated MW:

40kDa/54kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55223

Uniprot

Q9BVG3

Synonyms:

TRIM62; DEAR1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 120-280 of human TRIM62 (NP_060677.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

