

Product Information

Product SKU:	CAB15867	Gene ID:	55860	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in retrograde axonal transport of mitochondrion. Predicted to be located in cytosol; extracellular region; and secretory granule. Predicted to be part of dynactin complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	Arp10; Arp11; ACTR11; HARP11; ACTR10
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 218-417 of human ACTR10 (NP_060947.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.