

Product Information

Product SKU:	CAB15876	Gene ID:	57446	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	41kDa	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in signal transduction. Located in cytoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NDRG3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 184-363 of human NDRG3 (NP_071922.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.