

TRAF3IP3 Rabbit Polyclonal Antibody



CAB15899

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

40kDa/61kDa/63kDa

Applications:

WB

Reactivity:

Mouse, Rat

Protein Background

The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

80342

Uniprot

Q9Y228

Synonyms:

TRAF3IP3; T3JAM

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 290-530 of human TRAF3IP3 (NP_079504.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
