

Defb26 Rabbit Polyclonal Antibody



CAB15903

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

14kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

641629

Uniprot

Q32ZG6

Synonyms:

DEFB126; C20orf8; DEFB-26; DEFB26; HBD26; bA530N10.1; hBD-26


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 22-128 of rat Defb26 (NP_001032595.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of U2OS cells using Defb26 antibody (CAB15903) at dilution of 1:100. Blue: DAPI for nuclear staining.