

## CAB15949

---

### Product Information

<b>Product SKU:</b>	CAB15949	<b>Gene ID:</b>	127262	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	30kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	30kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable identical protein binding activity. Predicted to be involved in regulation of glutamatergic synaptic transmission. Located in extracellular exosome.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** TPRGL; mover; FAM79A; SVAP30; h-mover; TPRG1L

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human TPRG1L (NP\_877429.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.