

Product Information

Product SKU:	CAB15974	Gene ID:	284382	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be located in cytoplasm. Predicted to be part of dynactin complex. Implicated in spermatogenic failure 53.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	HSD21; SPGF53; ACTL9
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-416 of human ACTL9 (NP_848620.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.