

Product Information

Product SKU:	CAB15977	Gene ID:	337867	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	39kDa	Conjugate:	Unconjugated
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

**Background:** Involved in negative regulation of canonical Wnt signaling pathway and negative regulation of retrograde protein transport, ER to cytosol. Acts upstream of or within protein localization to endoplasmic reticulum. Located in endoplasmic reticulum.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PHGDHL1; UBAC2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human UBAC2 (NP_001137544.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.