

SRRM2 Rabbit Polyclonal Antibody



CAB15993

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

34kDa/256kDa/299kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23524

Uniprot

Q9UQ35

Synonyms:

SRRM2; 300-KD; CWF21; Cwc21; HSPC075; SRL300; SRm300

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRRM2 (NP_057417.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
