

C1QTNF7 Rabbit Polyclonal Antibody



CAB16031

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

30kDa/31kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

114905

Uniprot

Q9BXJ2

Synonyms:

C1QTNF7; CTRP7; ZACRP7

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 24-190 of human C1QTNF7 (NP_001128642.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
