

Product Information

Product SKU:	CAB16076	Gene ID:	8079	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in regulation of transcription, DNA-templated. Located in membrane.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:100 - 1:200
Synonyms:	NTN4; MLF2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human MLF2 (NP_005430.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.