

MLF2 Rabbit Polyclonal Antibody



CAB16076

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

28kDa

Calculated MW:

28kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC
1:100 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

8079

Uniprot

Q15773

Synonyms:

MLF2; NTN4


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human MLF2 (NP_005430.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images




Western blot analysis of extracts of various cell lines, using MLF2 antibody (CAB16076) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 5s.



Immunohistochemistry of paraffin-embedded rat liver using MLF2 antibody (CAB16076) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human esophageal using MLF2 antibody (CAB16076) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse spleen using MLF2 antibody (CAB16076) at dilution of 1:100 (40x lens).