

Product Information

Product SKU:	CAB16115	Gene ID:	27153	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	85kDa	Conjugate:	Unconjugated
Calculated MW:	94kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ZNF777
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 532-831 of human ZNF777 (NP_056509.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.