

# EBF4 Rabbit Polyclonal Antibody



CAB16139

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

64kDa/68kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

EBF4 belongs to the conserved Olf/EBF family of helix-loop-helix transcription factors, members of which play important roles in neural development and B-cell maturation.

## Immunogen information

### Gene ID:

57593

### Uniprot

Q9BQW3

### Synonyms:

EBF4; COE4; O/E-4

### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 50-150 of human EBF4 (NP\_001103984.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---