
Product Information

Product SKU:	CAB16145	Gene ID:	79838	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	114kDa	Conjugate:	Unconjugated
Calculated MW:	115kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable mechanosensitive ion channel activity. Predicted to be involved in ion transmembrane transport. Located in extracellular exosome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: TMC5

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 676-760 of human TMC5 (NP_079056.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.