

Product Information

Product SKU:	CAB16155	Gene ID:	84901	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in positive regulation of transcription by RNA polymerase II. Predicted to be located in cytoplasm and nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ESC2; NIP45; RAD60; NFATC2IP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 290-360 of human NFATC2IP (NP_116204.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.