

CAB16161

Product Information

Product SKU:	CAB16161	Gene ID:	116511	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable G protein-coupled receptor activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Located in cytosol; nucleoplasm; and plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MRG; MAS-L; dJ994E9.2; MAS1L
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human MAS1L (NP_443199.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.