

CAB16164

Product Information

Product SKU:	CAB16164	Gene ID:	121256	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	122kDa	Conjugate:	Unconjugated
Calculated MW:	122kDa	Isotype:	IgG

Immunogen Information

Background: May serve as a cell-surface marker for oligodendrocyte differentiation.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: TMEM132D; MOLT; PPP1R153

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 730-915 of human TMEM132D (NP_597705.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.