

CAB16192

Product Information

Product SKU:	CAB16192	Gene ID:	55110	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity. Involved in mRNA splicing, via spliceosome and nuclear-transcribed mRNA catabolic process, nonsense-mediated decay. Located in nucleus. Part of U2-type catalytic step 1 spliceosome; U2-type precatalytic spliceosome; and exon-exon junction complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MGN2; mago; magoh; MAGOHB
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-148 of human MAGOHB (NP_060518.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.