

CAB16215

Product Information

Product SKU:	CAB16215	Gene ID:	114990	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	71kDa	Conjugate:	Unconjugated
Calculated MW:	72kDa	Isotype:	IgG

Immunogen Information

Background: Enables transforming growth factor beta binding activity. Involved in negative regulation of epithelial to mesenchymal transition and negative regulation of transforming growth factor beta receptor signaling pathway. Located in cell surface and extracellular exosome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: SLITL2; VASN

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 420-575 of human VASN (NP_612449.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.