

NDUFA3 Rabbit Polyclonal Antibody



CAB16238

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

10kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4696

Uniprot

O95167

Synonyms:

NDUFA3; B9; CI-B9

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-84 of human NDUFA3 (NP_004533.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
