

Product Information

Product SKU:	CAB16330	Gene ID:	55323	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	54kDa	Conjugate:	Unconjugated
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity and myosin binding activity. Involved in positive regulation of collagen biosynthetic process; positive regulation of mRNA binding activity; and positive regulation of translation. Located in nucleus. Part of polysome.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ACHN; LARP6
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human LARP6 (NP_060827.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.