

ApoE Rabbit Polyclonal Antibody



CAB16344

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

36kDa

Calculated MW:

36kDa

Applications:

WB IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

11816

Uniprot

P02649

Synonyms:

APOE; AD2; APO-E; ApoE4; LDLCQ5; LPG; AI255918; Apo-E; Apoe

Immunogen:

A synthetic peptide of mouse Apoe

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

