

Product Information

Product SKU:	CAB16362	Gene ID:	20530	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	16kDa	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

Background:	Acts upstream of or within cellular copper ion homeostasis and regulation of copper ion transmembrane transport. Located in late endosome; membrane; and recycling endosome. Is expressed in brain; retina nuclear layer; and urinary system. Orthologous to human SLC31A2 (solute carrier family 31 member 2).
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CTR2; Slc31a2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse Slc31a2 (NP_001277447.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.