

# FOXA3 Rabbit Polyclonal Antibody



CAB16382

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

37kDa

### Applications:

WB

### Reactivity:

Rat

## Protein Background

This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. The crystal structure of a similar protein in rat has been resolved.

## Immunogen information

### Gene ID:

3171

### Uniprot

P55318

### Synonyms:

FOXA3; FKHH3; HNF3G; TCF3G

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXA3 (NP\_004488.2).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---