

PLAA Rabbit Polyclonal Antibody



CAB16466

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

87kDa

Applications:

IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9373

Uniprot

Q9Y263

Synonyms:

PLAA; DOA1; PLA2P; PLAP; NDMSBA


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 536-795 of human PLAA (NP_001026859.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded rat ovary using PLAA antibody (CAB16466) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human placenta using PLAA antibody (CAB16466) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using PLAA antibody (CAB16466) at dilution of 1:100 (40x lens).