

CAB16481

Product Information

Product SKU:	CAB16481	Gene ID:	10394	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa	Isotype:	IgG

Immunogen Information

**Background:** An extracellular matrix structural constituent conferring compression resistance. Involved in several processes, including granulocyte activation; histamine biosynthetic process; and regulation of gene expression. Located in collagen-containing extracellular matrix.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MBP2; MBPH; PRG3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 18-225 of human PRG3 (NP_006084.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.