

Product Information

Product SKU:	CAB16493	Gene ID:	11170	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to enable actin binding activity. Involved in several processes, including negative regulation of G1/S transition of mitotic cell cycle; negative regulation of focal adhesion assembly; and regulation of cytoskeleton organization. Located in several cellular components, including focal adhesion; ruffle membrane; and stress fiber.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	DRR1; TU3A; FAM107A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human FAM107A (NP_009108.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.