

CD2AP Rabbit Polyclonal Antibody



CAB16498

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

71kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

12488

Uniprot

Q9Y5K6

Synonyms:

Cd2ap; AL024079; C78928; METS-1; Mets1; CD2AP

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human CD2AP (NP_036252.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
