

BRPF3 Rabbit Polyclonal Antibody



CAB16515

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

99kDa/106kDa/135kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

27154

Uniprot

Q9ULD4

Synonyms:

BRPF3

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 50-150 of human BRPF3 (NP_056510.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
